



Modell Nr.	ML1309WN
Temperaturbereich	Kühlteil: +2 °C bis +10 °C Tiefkühlteil: -18 °C bis -25 °C
Umluft	nein
Abtauung (automatisch / manuell)	Kühlteil: automatisch Tiefkühlteil: manuell
Elektronikregler	ja
Digitalanzeige	ja
Temperaturalarm optisch / akustisch	ja
Über- und Untertemperaturalarm	ja
Türalarm	ja
Potentialfreier Alarmkontakt	ja
Außengehäuse	Stahlblech, lackiert
Farbe	weiß
Türschloss	ja
Türtyp	Volltür
Türanschlag	wechselbar
Türdichtung	wechselbar
Innenraum	Kunststoff, weiß
Innenbeleuchtung	nein
Nutzraumdurchführung Ø 16 mm	ja
„Innenraum frei von Zündquellen“	optional
Innenausstattung Kühlteil	
• Gitterauflagen (Anzahl)	1
• Auflagepositionen (Anzahl)	3
• Maße Gitterauflage (mm) (BxTxH)	495x325x6
• Schubladen (Anzahl)	0
• Drahtkörbe (Anzahl)	1
• Verstellbare Türablagen (Anzahl)	3
• Optionales Schubladensystem	ja
Innenausstattung Tiefkühlteil	
• Schubladen (Anzahl)	0
Verstellbare Füße (Anzahl & Position)	4 Füße 2 Rollen hinten
Inhalt (l), brutto	138 120 Kühlteil / 18 Tiefkühlteil
Inhalt (l), netto	136 118 Kühlteil / 18 Tiefkühlteil
Breite (mm)	595
Tiefe (mm)	600
Tiefe bei geöffneter Tür (mm)	1150
Höhe (mm)	850
Innenmaße (mm) (BxTxH)	Kühlteil: 500x480x480 Tiefkühlteil: 445x290x160 Stufe (TxH): 155x215
Gewicht (kg)	34
Netzanschluss	1/N/PE-230 V 50 Hz
Versandbreite (mm)	700
Versandtiefe (mm)	750
Versandhöhe (mm)	1027
Versandgewicht (kg)	40
Kältemittel FCKW- / HFCKW-frei	R600a



ML1309WN

Passende Ordnungssystemvorschläge finden Sie in unserem separaten Katalog. Bitte fordern Sie Ihren persönlichen Einrichtungsvorschlag an.





Model No.	ML1309WN
Temperature range	Fridge: +2 °C to +10 °C Freezer: -18 °C to -25 °C
Forced air circulation	no
Defrost (automatic / manual)	Fridge: automatic Freezer: manual
Electronic controller	yes
Digital display	yes
Temperature alarm visual / audible	yes
Excess and insufficient temp. alarm	yes
Open door alarm	yes
No-volt remote alarm	yes
External finish	painted steel plate
Colour	white
Lock	yes
Door type	solid door
Door hinges	changeable
Door seal	changeable
Internal finish	plastic, white
Interior light	no
Port hole Ø 16 mm	yes
Spark free	optional
Interior fittings Fridge	
• Shelves (number)	1
• Shelf positions (number)	3
• Shelf size (mm) (WxDxH)	495x325x6
• Drawers (number)	0
• Wire baskets (number)	1
• Adjustable door bins (number)	3
• Optional drawer system	Yes
Interior fittings Freezer	
• Drawers (number)	0
Adjustable feet (number & position)	4 feet 2 rear rollers
Gross capacity (l)	138 120 Fridge / 18 Freezer
Net capacity (l)	136 118 Fridge / 18 Freezer
External width (mm)	595
External depth (mm)	600
Door open depth (mm)	1150
External height (mm)	850
Internal dimensions (mm) (DxWxH)	Fridge: 500x480x480 Freezer: 445x290x160 step (TxH): 155x215
Product weight (kg)	34
Voltage supply	1/NPE-230 V 50 Hz
Packed width (mm)	700
Packed depth (mm)	750
Packed height (mm)	1027
Packed weight (kg)	40
Refrigerant type CFC / HCFC free	R600a



ML1309WN

Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

