



Model No.	ML0521WN	ML1225WN	ML1625WN	ML3125WN
Temperature range	-18°C to -25°C	-18°C to -25°C	-18°C to -25°C	-18°C to -25°C
Forced air circulation	No	No	No	No
Defrost (automatic / manual)	manual	manual	manual	manual
Electronic controller	Yes	Yes	Yes	Yes
Digital display	Yes	Yes	Yes	Yes
Temperature alarm visual / audible	Yes	Yes	Yes	Yes
Excess and insufficient temp. alarm	Yes	Yes	Yes	Yes
Open door alarm	No	No	Yes	Yes
No-volt remote alarm	No	No	Yes	Yes
External finish	painted steel plate	painted steel plate	painted steel plate	painted steel plate
Colour	white	white	white	white
Lock	Yes	Yes	Yes	Yes
Door type	solid door	solid door	solid door	solid door
Door hinges	changeable	changeable	changeable	changeable
Door seal	not changeable	not changeable	not changeable	not changeable
Internal finish	plastic, white	plastic, white	plastic, white	plastic, white
Interior light	No	No	No	No
Spark free	Yes	Yes	Yes	Yes
Shelves (number)	1	0	0	0
Cooling spirals (number)	0	2	2	5
Shelf positions (number)	1	2 (fixed)	2 (fixed)	5 (fixed)
Shelf size (mm) (WxDxH)	340x180x6	-	-	-
Drawers (number)	0	0	0	4 and 2 flaps
Wire baskets (number)	0	1	3	0
Optional drawer system	Yes	Yes	Yes	Yes
Adjustable feet (number & position)	No	1 front foot	2 front feet 2 rear rollers	2 front feet 2 rear rollers
Gross capacity (l)	50	102	151	232
Net capacity (l)	38	87	134	212
External width (mm)	500	495	595	595
External depth (mm)	510	620	595	634
Door open depth (mm)	976	1050	1165	1185
External height (mm)	562	830	845	1570
Internal dimensions (mm) (WxDxH)	345x330x365	414x415x683	500x450x705 step (DxH): 165x185	432x430x1300 step (DxH): 190x255
Product weight (kg)	17	33	36	60
Voltage supply	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz	1/N/PE-230V 50Hz
Packed width (mm)	700	700	700	700
Packed depth (mm)	750	750	750	750
Packed height (mm)	740	1000	1022	1735
Packed weight (kg)	24	40	43	67
Refrigerant type CFC / HCFC free	R134a	R134a	R600a	R600a



ML0521WN



ML0521WN with Inventory System (optional)



ML1225WN



ML1225WN with Inventory System (optional)



ML1625WN



ML3125WN



Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

