



Model No.	ML1307WN	ML2726WN	ML3226WN
Temperature range	Fridge: +2 °C to +10 °C Freezer: -18 °C to -20 °C	Fridge: +2 °C to +10 °C Freezer: -18 °C to -20 °C	Fridge: +2 °C to +10 °C Freezer: -18 °C to -20 °C
Forced air circulation	No	No	No
Defrost (automatic / manual)	Fridge: automatic Freezer: manual	Fridge: automatic Freezer: manual	Fridge: automatic Freezer: manual
Electronic controller	Yes	Yes	Yes
Digital display	Yes	Yes	Yes
Temperature alarm visual / audible	Yes	Yes	Yes
Excess and insufficient temp. alarm	Yes	Yes	Yes
Open door alarm	Yes	Yes	Yes
No-volt remote alarm	Yes	Yes	Yes
External finish	painted steel plate	painted steel plate	painted steel plate
Colour	white	white	white
Lock	Yes	Yes (2)	Yes (2)
Door type	solid door	solid door	solid door
Door hinges	changeable	changeable	changeable
Door seal	not changeable	not changeable	not changeable
Internal finish	plastic, white	plastic, white	plastic, white
Interior light	No	No	No
Spark free	optional	optional	optional
Interior fittings Fridge			
• Shelves (number)	2	3	3
• Shelf positions (number)	4	6	6
• Shelf size (mm) (WxDxH)	1: 448x308x4.5 1: 448x200x4.5	448x308x4.5	505x320x6
• Drawers (number)	1	2 half	2 half
• Wire baskets (number)	0	0	0
• Adjustable door bins (number)	2 full and 2 half	2 full and 2 half	2 full and 4 half
• Optional drawer system	Yes	Yes	Yes
Interior fittings Freezer			
• Drawers (number)	0	2 and 1 flap	2 and 2 flaps
Adjustable feet (number & position)	2 front feet 2 rear rollers	2 front feet 2 rear rollers	2 front feet 2 rear rollers
Gross capacity (l)	119 106 Fridge / 13 Freezer	263 159 Fridge / 104 Freezer	324 190 Fridge / 134 Freezer
Net capacity (l)	117 104 Fridge / 13 Freezer	245 148 Fridge / 97 Freezer	302 177 Fridge / 125 Freezer
External width (mm)	540	540	600
External depth (mm)	600	600	600
Door open depth (mm)	1090	1110	1190
External height (mm)	860	1680	1870
Internal dimensions (mm) (DxWxH)	Fridge: 450x450x525 Freezer: 400x295x120 step (TxH): 155x225	Fridge: 457x450x725 Freezer: 408x415x670 step (TxH): 170x220	Fridge: 510x430x805 Freezer: 463x445x780 step (TxH): 250x260
Product weight (kg)	34	58	70
Voltage supply	1/N/PE-230 V 50 Hz	1/N/PE-230 V 50 Hz	1/N/PE-230 V 50 Hz
Packed width (mm)	700	700	700
Packed depth (mm)	750	750	750
Packed height (mm)	1027	1850	2035
Packed weight (kg)	41	65	77
Refrigerant type CFC / HCFC free	R600a	R600a	R600a



ML1307WN



ML2726WN

ML3226WN

Recommendations for additional storage sets can be found in our separate catalogue. Please ask for a personal fitting option.

