

Small Scale Chest Freezer ProfiMaster *Mini*

- Static cooling and manual defrost
- Low-noise refrigeration system, air-cooled
- Operating panel with electronic control and digital display (LED)
 - Visual/audible alarm device
 - Alarm warning for excess and insufficient temperature
- External housing made of stainless steel (CNS)
- Interior made of stainless steel (CNS)
- Insulated lid with hinge and magnetic seal



Small Scale Chest Freezer PMM0860EEN

Model ProfiMaster **PMM0860EEN**

Temperature range	°C: -20/-50
	: Other temperature ranges available on request.
Content	litres (gross): 8
Ext. dimensions (WxDxH)	mm: 370x500x580
Int. dimensions (WxDxH)	mm: 150x300x200
Weight	kgs: 35
Mains connection	: 1/N/PE ~ 230 V 50 Hz
Refrigerant CFC- and HCFC-free:	R507

ProfiMaster Mini + Accessories **Outside / Inside** **Cat. No.**

ProfiMaster Mini PMM0860	CNS / CNS	GTPMM0860
No-volt remote alarm contact for connection to a remote alarm system.		ZPMM2-009
ProfiTec Temperature Monitoring and Recording System (see sep. brochure).		GPTEC4A
ProfiCaller telephone dialler (see sep. brochure).		on request
Pt-100 temperature sensor, add. installed for ext. temperature recording. Forced air circulation.		ZPM2-023-2L on request
MIN/MAX thermometer to be placed in refrigerators or freezers. Measuring range -50/+50 °C, an add. temperature sensor is not required.		ZXX2-019-001
MIN/MAX digital thermometer. Hand-held unit with digital display as well as data-hold and relative measuring function. Measuring range -200/+1370 °C. An add. temperature sensor is required, eg. Cat.-No. ZXX2-022-100.		ZXX2-019-011
Temperature sensor for MIN/MAX thermometer ZXX2-019-011. Measuring range -100/+1000 °C. Length 2 m. Other ranges on request.		ZXX2-022-100

Innovation prize

ProfiTec Intelligent Temperature Monitoring and Recording System



- For up to 4 refrigerators and/or freezers
- Readings are stored over a period up to 781 days
- Warning for excess and insufficient temperature
- Includes No-volt remote alarm contact
- The measuring protocol can be fed to PC via RS 232 port



MIN/MAX digital thermometer with add. temperature sensor

Laboratory Upright Freezer ProfiLine Pegasus PLPE0186EWN



Temperature range	°C: -40/-90
Content	litres (gross): 7
Ext. dimensions (WxDxH)	mm: 400x630x665
Int. dimensions (WxDxH)	mm: 150x143x310
Depth with door opened	mm: 960
Weight	kgs: 54
Mains connection	: 1/N/PE ~ 230 V 50-60 Hz

See separate brochure.